




Recombinant Human Peroxiredoxin-like 2A (PRXL2A)

Product Code	CSB-EP883589HU
Abbreviation	Recombinant Human PRXL2A protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BRX8
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSFLQDPSFFTGMGMWSIGAGALGAAALALLLANTDVFLSKPQKAALEYLEDIDL KTLEKEPRTFKAKELWEKNGAVIMAVRRPGCFLCREEAADLSSLKSMLDQLGV PLYAVVKEHIRTEVKDFQPYFKGEIFLDEKKKFYGPQRRKMMFMGFIRLGVMWY NFFRAWNGGFSGNLEGEFILGGVFVVGSGKQGILLEHREKEFGDKVNLLSVL EAAKMIKPQTLASEKK
Research Area	Cancer?
Source	E.coli
Target Names	FAM213A
Protein Names	Recommended name: Redox-regulatory protein FAM213AAlternative name(s): Peroxiredoxin-like 2 activated in M-CSF stimulated monocytes Short name= Protein PAMM
Expression Region	1-229aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	33.2 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.