




# Recombinant Mouse Thyrotropin receptor (Tshr), partial

<b>Product Code</b>	CSB-BP025131MO
<b>Abbreviation</b>	Recombinant Mouse Tshr protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P47750
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	KECASPPCECHQEDDFRVTCKELHRIPSLPPSTQTLKLIETHLKTIPSLAFSSLP NISRIYLSIDATLQRLEPHSFYNLSKMTHEIRNTRSLTYIDPDALTEPLLKFLGIF NTGLRIFPDLTKIYSTDIFFILEITDNPYMTSVPENAFQGLCNETLTLKLYNNGFT SVQGHAFNGTKLDAVYLNKNKYLT AIDNDAFGGVYSGPTLLDVSSTSVTALPS KGLEHLKELIAKDTWTLKKLPLSLSFLHLTRADLSYPHCCAFKNQKKIRGILES LMCNESSIRNLRQRKSVNILRGPIYQEYEEDPGDNSVGKQNSKFQESPSNSH YYVFFEEQEDEVVGFQELKNPQEETLQAFESHYDYTVCGDNEDMVCTPKSD EFNPCEDIMG
<b>Research Area</b>	Others
<b>Source</b>	Baculovirus
<b>Target Names</b>	Tshr
<b>Protein Names</b>	Recommended name: Thyrotropin receptor Alternative name(s): Thyroid-stimulating hormone receptor Short name= TSH-R
<b>Expression Region</b>	22-413aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	45.6
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.