




Recombinant Mouse Thyrotropin receptor (Tshr), partial

Product Code	CSB-MP025131MO
Abbreviation	Recombinant Mouse Tshr protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47750
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	KECASPPCECHQEDDFRVTCKELHRIPSLPPSTQTLKLIETHLKTIPSLAFSSLP NISRIYLSIDATLQRLEPHSFYNLSKMTHEIRNTRSLTYIDPDALTELPLLKFLGIF NTGLRIFPDLTKIYSTDIFFILEITDNPYMTSVPENAFQGLCNETLTLKLYNNGFT SVQGHAFNGTKLDAVYLNKNKYLT AIDNDAFGGVYSGPTLLDVSSTSVTALPS KGLEHLKELIAKDTWTLKKLPLSLSFLHLTRADLSYP SHCCAFKNQKKIRGILES LMCNESSIRNLRQRKSVNLRGPIYQEYEEDPGDNSVG YKQNSKFQESPSNSH YYVFFEEQEDEVVGFGQELKNPQEETLQAFESHYDYTVCGDNEDMVCTPKSD EFNPCEDIMG
Research Area	Others
Source	Mammalian cell
Target Names	Tshr
Protein Names	Recommended name: Thyrotropin receptor Alternative name(s): Thyroid-stimulating hormone receptor Short name= TSH-R
Expression Region	22-413aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	46.9
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.