




Recombinant Human Probable histone-lysine N-methyltransferase PRDM7 (PRDM7)

Product Code	CSB-EP873626HU
Abbreviation	Recombinant Human PRDM7 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NQW5
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSPERSQEE SPEGDTERTERKPMVKDAFKDISIYFTKEEWAEMGDWEKTRYR NVKMNYNALITVGLRATRP AFMCHRRQA IKLQVDDTEDSDEEWTPRQQVKPP WMAFRGEQSKHQKGMPKASFNNESLRELSGTPNLLNTSDSEQAQKPVSP GEASTSGQHSRLKLELRRKETEGKMYSLRERKGHAYKEISEPQDDDYLYCEM CQNFIDSCAAHGPPTFVKDSAVDKGHPNRSALSLPPGLRIGPSGIPQAGLG VWNEASDLPLGLHFGPYEGRITDEEAAANS GYSWLITKGRNCYEYVDGKDKSS ANWMRYVNCARDDEEQNLVAFQYHRQIF YRTCRVIRPGCELLVWSGDEYGG ELGIRSSIEPAESLGQAVNCWSGMGMSMARNWASSGAASGRKSSWQGENQ SQRSIHVPHAVWPFQVKNFSVNMWNAITPLRTSQDHLQENFSNQRI PAQGIRI RSGNILIHA AVMTKPKVKRSKKGPN S
Research Area	Others
Source	E.coli
Target Names	PRDM7
Protein Names	Recommended name: Probable histone-lysine N-methyltransferase PRDM7 EC= 2.1.1.43 Alternative name(s): PR domain zinc finger protein 7 PR domain- containing protein 7
Expression Region	1-492aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	59.9 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.