




# Recombinant Human Heterogeneous nuclear ribonucleoprotein H2 (HNRNPH2)

<b>Product Code</b>	CSB-EP346289HU
<b>Abbreviation</b>	Recombinant Human HNRNPH2 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P55795
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	MMLSTEGREGFVVKVRGLPWSCSADEVMRFFSDCKIQNGTSGIRFIYTREGR PSGEAFVELESEEEVVKLALKKDRETMGHRYVEVFKSNSVEMDWVLKHTGPNS PDTANDGFVRLRGLPFGCSKEEIVQFFSGLIIVPNGMTLPVDFQGRSTGEAFV QFASQEIAEKALKKKHKERIGHRYIEIFKSSRAEVRTHYDPPRKLMMAMQRPQPYD RPGAGRGYNSIGRGAGFERMRRGAYGGGYGGYDDYGGYNDGYGFGSDRF GRDLNYCFSGMSDHRYGDGGSSFQSTTGHCVHMRGLPYRATENDIYNFFSP LNPMRVHIEIGPDGRVTGEADVEFATHEDAVAAMAKDKANMQHRYVELFLNS TAGTSGGAYDHSYVELFLNSTAGASGGAYGSQMMGGMGLSNQSSYGGPAS QQLSGGYGGGYGGQSSMSGYDQVLQENSDDYQSNLA
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	HNRNPH2
<b>Protein Names</b>	Recommended name: Heterogeneous nuclear ribonucleoprotein H2 Short name= hnRNP H2 Alternative name(s): FTP-3 Heterogeneous nuclear ribonucleoprotein H' Short name= hnRNP H'
<b>Expression Region</b>	1-449aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	56.7 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.