




Recombinant Rat Nuclear pore glycoprotein p62 (Nup62)

Product Code	CSB-EP016204RAb0
Abbreviation	Recombinant Rat Nup62 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17955
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	SGFNFGGTGAPAGGFTFGTAKTATTTTPATGFSFSASGTGTGGFNFGTPSQPA ATTPSTSLFSLATQTSTTQTPGFNFGTTPASGGTGFSLGISTPKLSLSSTAATP ATANTGSFGLGSSTLTNAISGASTSSQGTAPTGFVFGSSTTSAPSTGTTGFSF TSGSASQPGASGFNIGSVGSLAQPTALSGSPFTPATLATTTAGATQPAAATPT AATTSAGSTLFASIAAAPASSSTTVLSLSAPATTAATPTAGTLGFSLKAPGAAP GASTTSTTTTTTTTTTTASTSSSTTTTGFAISLKPLVPAGPSSVAATALPASSTA VGTTTGPMATYAQLESINKWSLELEDQERHFLQQATQVNAWDRTLIENGEKI TSLHREVEKVKLDQKRLDQELDFILSQQKELEDLLSPLEESVKEQSGTIYQLQA DEEREKTYKLAENIDAQLKRMAQDLKDIIHLNMAGGPADTSDPLQQICKILNA HMDSLQWVDQSSALLQRRVEEASRVCESRRKEQERSLRIAFD
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	Nup62
Expression Region	2-525aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	59.3 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.