




Recombinant Human Fibrinogen-like protein 1 (FGL1)

Product Code	CSB-MP008653HU
Abbreviation	Recombinant Human FGL1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q08830
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 95% as determined by SDS-PAGE
Sequence	LEDCAQEQMRLRAQVRLLETRVKQQQVKIKQLLQENEVQFLDKGDENTVIDLG SKRQYADCSEIFNDGYKLSGFYKIKPLQSPAEFVSVYCDMSDGGGWTVIQRRS DGSENFNRGWKDYENGFVQKHGEYWLGNKNLHFLTTQEDYTLKIDLADF EKNSRYAQYKNFKVGDEKNFYELNIGEYSGTAGDSLGNFHPFVQWWASHQ RMKFSTWDRDHDNYEGNCAEEDQSGWVFNRRCHSANLNGVYYSGPYTAKTD NGIVWYTWHGWWYSLKSVVMKIRPNDPFI
Research Area	Cardiovascular?
Source	Mammalian cell
Target Names	FGL1
Expression Region	23-312aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal hFc1-Myc-tagged
Mol. Weight	64.1 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.