




Recombinant Human Tyrosine-protein phosphatase non-receptor type substrate 1 (SIRPA), partial

Product Code	CSB-MP021334HUh6
Abbreviation	Recombinant Human SIRPA protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P78324
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	EEELQVIQPDKSVLVAAGETATLRCTATSLIPVGPIQWFRGAGPGRELIYNQKE GHFPRVTTVSDLTKRNNMDFSIRIGNITPADAGTYCYVKFRKGSPDDVEFKSG AGTELSVRAKPSAPVWSGPAARATPQHTVSFTCESHGFSPRDITLKWFKNGN ELSDFQTNVDPVGESVSYSIHSTAKVVLTREDVHSQVICEVAHVTLQGDPLRG TANLSETIRVPPTLEVTQQPVRAENQVNVTCQVRKFYPQRLQLTWLENGNVS RTETASTVTENKDGTYNWMSWLLNVNSAHRDDVKLTCQVEHDGQPAVSKSH DLKVS AHPKEQGSNTAAENTGSNER
Research Area	Cardiovascular
Source	Mammalian cell
Target Names	SIRPA
Expression Region	31-370aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal hFc1-Myc-tagged
Mol. Weight	67.40000000000001
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.