




Recombinant Human 40S ribosomal protein S9 (RPS9)

Product Code	CSB-BP020482HU
Abbreviation	Recombinant Human RPS9 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46781
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	PVARSWVCRKTYVTPRRPF EKSR LDQELKLIGEYGLRNKREVWRVKFTLAKIR KAARELLTLDEKDP RRLFEGNALLRRLVRIGVLDEGKMKLDYILGLKIEDFLERR LQTQVFKLGLAKSIHHARVLIRQRHIRVRKQVNVNIPSFIVRLDSQKHIDFSLRSPY GGGRPGRVKRKNNAKKGQGGAGAGDDEEED
Research Area	Epigenetics and Nuclear Signaling
Source	Baculovirus
Target Names	RPS9
Protein Names	Recommended name: 40S ribosomal protein S9
Expression Region	2-194aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	26.3
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.