




# Recombinant Leiurus hebraeus Alpha-like toxin Lqh6

|                          |                                                                                                                                                                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP348580LDS                                                                                                                                                                                                                                                                 |
| <b>Abbreviation</b>      | Recombinant Leiurus hebraeus Alpha-like toxin Lqh6 protein                                                                                                                                                                                                                      |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| <b>Uniprot No.</b>       | P59356                                                                                                                                                                                                                                                                          |
| <b>Product Type</b>      | Recombinant Proteins                                                                                                                                                                                                                                                            |
| <b>Immunogen Species</b> | Leiurus hebraeus (Deathstalker scorpion) (Leiurus quinquestriatus hebraeus)                                                                                                                                                                                                     |
| <b>Purity</b>            | ≥ 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                                |
| <b>Sequence</b>          | VRDGYIAQPENCYVYHCIPDCDTLCKDNGGTGGHCGFKLGHGIACWCNALPDN<br>VGIIVDGVKCHK                                                                                                                                                                                                           |
| <b>Research Area</b>     | Others                                                                                                                                                                                                                                                                          |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                     |
| <b>Target Names</b>      | N/A                                                                                                                                                                                                                                                                             |
| <b>Protein Names</b>     | Recommended name: Alpha-like toxin Lqh6<br>Alternative name(s): Lqh VI<br>Short name= LqhVI                                                                                                                                                                                     |
| <b>Expression Region</b> | 1-64aa                                                                                                                                                                                                                                                                          |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                             |
| <b>Tag Info</b>          | N-terminal 10xHis-tagged and C-terminal Myc-tagged                                                                                                                                                                                                                              |
| <b>Mol. Weight</b>       | 10.7                                                                                                                                                                                                                                                                            |
| <b>Protein Length</b>    | Full Length                                                                                                                                                                                                                                                                     |

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.