




Recombinant Human Chymotrypsin-like elastase family member 2B (CELA2B)

Product Code	CSB-BP362355HU
Abbreviation	Recombinant Human CELA2B protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08218
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MLGGEEARPNSWPWQVSLQYSSNGQWYHTCGGSLIANSWVLTAAHCISSSG IYRVMLGQHNLVVAESGLAVSVSKIVVHKDWNSDQVSKGNDIALLLKLANPVS LTDKIQLACLPPAGTILPNNYPCYVTGWGRLQTNGALPDDLKQGQLLVVDYAT CSSSGWWGSTVKTNMICAGGDGVICTCNGDSGGPLNCQASDGRWEVHGIG SLTSVLGCNYYYKPSIFTRVSNYNDWINSVIANN
Research Area	Cell Biology
Source	Baculovirus
Target Names	CELA2B
Protein Names	Recommended name: Chymotrypsin-like elastase family member 2B EC= 3.4.21.71 Alternative name(s): Elastase-2B
Expression Region	29-269aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	28.4
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.