




Recombinant Human Arylacetamide deacetylase (AADAC), partial

Product Code	CSB-EP001014HUb0
Abbreviation	Recombinant Human AADAC protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P22760
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	PDNVEEPWRMMWINAHLKTIQNLATFVELLGLHHFMDSFKVVGSFDEVPTS DENVTVTETKFNNILVRVYVPKRKSEALRRGLFYIHGGGWCVGSAALSGYDLL SRWTADRLDAVVVSTNYRLAPKYHFPIQFEDVYNALRWFLRKKVLAKYGVNPE RIGISGDSAGGNLAAAVTQQLLDDPDVKIKLKIQSLIYPALQPLDLDLPSYQENS NFLFLSKSLMVRFWSEYFTTDRSLEKAMLSRQHVPVESHLLFKFVNWSSLLPE RFIKGHVYNNPNYGSSELAKKYPGFLDVRAAPLLADDNKLRLPLTYVITCQYD LLRDDGLMYVTRLRNTGVQVTHNHVEDGFHGAFSFLGLKISHRLINQYIEWLK ENL
Research Area	Signal Transduction
Source	E.coli
Target Names	AADAC
Expression Region	24-399aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	49.1 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.