




Recombinant Cricetulus longicaudatus Glutathione S-transferase P (GSTP1)

Product Code	CSB-EP009989DXV
Abbreviation	Recombinant Cricetulus longicaudatus GSTP1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46424
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Cricetulus longicaudatus(Long-tailed dwarf hamster)(Chinese hamster)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MPPYTIVYFPVRGRCEAMRILLADQGQSWKEEVVTVETWRKGSCLKSTCLYGQ LPKFEDGDLTYQSNAILRHLGRSLGLYGKDKQREAAALVDMVNDGVEDLRCKYI TLIYTKYEEGKDDYVKALPGHLKPFETLLSQNQGGKAFIVGDQISFVDYNLLDLL LIHQVLAPGCLDNFPLLSAYVARLSARPKIKAFLLSSPDHVNRPINGNGKQ
Research Area	Metabolism
Source	E.coli
Target Names	GSTP1
Protein Names	Recommended name: Glutathione S-transferase P EC= 2.5.1.18Alternative name(s): GST class-pi
Expression Region	1-210aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	28.7 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.