



# DPEP1 Antibody, Biotin conjugated

Product Code	CSB-PA007123LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16444
Immunogen	Recombinant Human Dipeptidase 1 protein (23-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dipeptidase 1 (EC 3.4.13.19) (Dehydropeptidase-I) (Microsomal dipeptidase) (Renal dipeptidase) (hRDP), DPEP1, MDP RDP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DPEP1