



DPEP1 Antibody, Biotin conjugated

Product Code	CSB-PA007123LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16444
Immunogen	Recombinant Human Dipeptidase 1 protein (23-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dipeptidase 1 (EC 3.4.13.19) (Dehydropeptidase-I) (Microsomal dipeptidase) (Renal dipeptidase) (hRDP), DPEP1, MDP RDP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DPEP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.