



INSRR Antibody, FITC conjugated

Product Code	CSB-PA011755LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14616
Immunogen	Recombinant Human Insulin receptor-related protein (747-921AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Insulin receptor-related protein (IRR) (EC 2.7.10.1) (IR-related receptor) [Cleaved into: Insulin receptor-related protein alpha chain, Insulin receptor-related protein beta chain], INSRR, IRR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INSRR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.