




EXOC3 Antibody

Product Code	CSB-PA527663LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60645
Immunogen	Recombinant Human Exocyst complex component 3 protein (410-756AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exocyst complex component 3 (Exocyst complex component Sec6), EXOC3, SEC6 SEC6L1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	EXOC3

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, MCF-7 whole cell lysate, A549 whole cell lysate, HL60 whole cell lysate, PC-3 whole cell lysate, U87 whole cell lysate

All lanes: EXOC3 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 87, 74, 86 kDa

Observed band size: 87 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.