



EXOC3 Antibody, HRP conjugated

Product Code	CSB-PA527663LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60645
Immunogen	Recombinant Human Exocyst complex component 3 protein (410-756AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exocyst complex component 3 (Exocyst complex component Sec6), EXOC3, SEC6 SEC6L1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	EXOC3