



EXOC3 Antibody, Biotin conjugated

Product Code	CSB-PA527663LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60645
Immunogen	Recombinant Human Exocyst complex component 3 protein (410-756AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exocyst complex component 3 (Exocyst complex component Sec6), EXOC3, SEC6 SEC6L1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	EXOC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.