



KCNMB3 Antibody, HRP conjugated

Product Code	CSB-PA878838LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPA1
Immunogen	Recombinant Human Calcium-activated potassium channel subunit beta-3 protein (82-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-activated potassium channel subunit beta-3 (BK channel subunit beta-3) (BKbeta3) (Hbeta3) (Calcium-activated potassium channel, subfamily M subunit beta-3) (Charybdotoxin receptor subunit beta-3) (K(VCA)beta-3) (Maxi K channel subunit beta-3) (Slo-beta-3), KCNMB3, KCNMB2 KCNMBL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KCNMB3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.