



Copper-metallothionein Antibody, Biotin conjugated

Product Code	CSB-PA345921ED01HDB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55947
Immunogen	Recombinant Helix pomatia Copper-metallothionein protein (1-64aa)
Raised In	Copper-metallothionein
Species Reactivity	Roman snail
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Copper-metallothionein?Cu-MT
Immunogen Species	Helix pomatia (Roman snail) (Edible snail)
Research Area	Others
Target Names	Copper-metallothionein
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.