



TJP1 Antibody, FITC conjugated

Product Code	CSB-PA023573HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07157
Immunogen	Recombinant Human Tight junction protein ZO-1 protein (116-414aa)
Raised In	Tight junction protein ZO-1
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tight junction protein ZO-1 (Tight junction protein 1) (Zona occludens protein 1) (Zonula occludens protein 1), TJP1, ZO1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TJP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.