



SEZ6L2 Antibody, Biotin conjugated

Product Code	CSB-PA761499LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXD5
Immunogen	Recombinant Human Seizure 6-like protein 2 protein (28-372aa)
Raised In	Seizure 6-like protein 2
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Seizure 6-like protein 2, SEZ6L2, PSK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SEZ6L2