




Fgf2 Antibody

Product Code	CSB-PA319495LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13109
Immunogen	Recombinant Rat Fibroblast growth factor 2 protein (10-154AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 2 (FGF-2) (Basic fibroblast growth factor) (bFGF) (Heparin-binding growth factor 2) (HBGF-2), Fgf2, Fgf-2
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Fgf2

Image



Western Blot

Positive WB detected in: Rat Heart tissue, Rat Brain tissue

All lanes: Fgf2 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 18 kDa

Observed band size: 18 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.