



# Fgf2 Antibody?FITC conjugated

<b>Product Code</b>	CSB-PA319495LC01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P13109
<b>Immunogen</b>	Recombinant Rat Fibroblast growth factor 2 protein (10-154AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Antigen Affinity purified & Affinity purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Fibroblast growth factor 2 (FGF-2) (Basic fibroblast growth factor) (bFGF) (Heparin-binding growth factor 2) (HBGF-2), Fgf2, Fgf-2
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	Fgf2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.