



IKBIP Antibody, HRP conjugated

Product Code	CSB-PA751229LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70UQ0
Immunogen	Recombinant Human Plexin-B2 protein (121-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inhibitor of nuclear factor kappa-B kinase-interacting protein (I kappa-B kinase-interacting protein) (IKBKB-interacting protein) (IKK-interacting protein), IKBIP, IKIP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IKBIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.