




# ZNF346 Antibody

<b>Product Code</b>	CSB-PA026691EA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UL40
<b>Immunogen</b>	Recombinant Human Angiopoietin-related protein 4 protein (1-294AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Antigen Affinity purified & Affinity purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Zinc finger protein 346 (Just another zinc finger protein), ZNF346, JAZ
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	ZNF346

## Image



### Western Blot

Positive WB detected in: HeLa whole cell lysate, A549 whole cell lysate, MCF-7 whole cell lysate, 293T whole cell lysate, Jurkat whole cell lysate, PC-3 whole cell lysate  
 All lanes: ZNF346 antibody at 1:1000  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/50000 dilution  
 Predicted band size: 33, 36, 30 kDa  
 Observed band size: 33 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.