



ZNF346 Antibody?FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA026691EC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UL40 |
| Immunogen | Recombinant Human Angiopoietin-related protein 4 protein (1-294AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Antigen Affinity purified & Affinity purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Zinc finger protein 346 (Just another zinc finger protein), ZNF346, JAZ |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | ZNF346 |