



ZNF346 Antibody? Biotin conjugated

Product Code	CSB-PA026691ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL40
Immunogen	Recombinant Human Angiopoietin-related protein 4 protein (1-294AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 346 (Just another zinc finger protein), ZNF346, JAZ
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF346
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.