




ANGPTL4 Antibody

Product Code	CSB-PA866314EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY76
Immunogen	Recombinant Human Trefoil factor 3 protein (26-406AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiopoietin-related protein 4 (Angiopoietin-like protein 4) (Hepatic fibrinogen/angiopoietin-related protein) (HFARP), ANGPTL4, ARP4 HFARP PGAR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ANGPTL4

Image



Western Blot

Positive WB detected in: SH-SY5Y whole cell lysate, HeLa whole cell lysate, A549 whole cell lysate

All lanes: ANGPTL4 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 46, 41, 27 kDa

Observed band size: 46 kDa