



ANGPTL4 Antibody, HRP conjugated

Product Code	CSB-PA866314EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY76
Immunogen	Recombinant Human Trefoil factor 3 protein (26-406AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiopoietin-related protein 4 (Angiopoietin-like protein 4) (Hepatic fibrinogen/angiopoietin-related protein) (HFARP), ANGPTL4, ARP4 HFARP PGAR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ANGPTL4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.