



# Mif Antibody? Biotin conjugated

|                            |                                                                                                                                                                                                                                                 |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA013826LD01MO                                                                                                                                                                                                                              |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                   |
| <b>Uniprot No.</b>         | P34884                                                                                                                                                                                                                                          |
| <b>Immunogen</b>           | Recombinant Mouse Macrophage migration inhibitory factor protein (2-115AA)                                                                                                                                                                      |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                                                                          |
| <b>Species Reactivity</b>  | Mouse                                                                                                                                                                                                                                           |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                                           |
| <b>Form</b>                | Liquid                                                                                                                                                                                                                                          |
| <b>Conjugate</b>           | Biotin                                                                                                                                                                                                                                          |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                                                                                                                                                                |
| <b>Purification Method</b> | >95%, Antigen Affinity purified & Affinity purified                                                                                                                                                                                             |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                      |
| <b>Alias</b>               | Macrophage migration inhibitory factor (MIF) (EC 5.3.2.1) (Delayed early response protein 6) (DER6) (Glycosylation-inhibiting factor) (GIF) (L-dopachrome isomerase) (L-dopachrome tautomerase) (EC 5.3.3.12) (Phenylpyruvate tautomerase), Mif |
| <b>Immunogen Species</b>   | Mus musculus (Mouse)                                                                                                                                                                                                                            |
| <b>Research Area</b>       | Cardiovascular                                                                                                                                                                                                                                  |
| <b>Target Names</b>        | Mif                                                                                                                                                                                                                                             |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                                                                                                                                                     |