



IFRD1 Antibody?FITC conjugated

Product Code	CSB-PA011069LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00458
Immunogen	Recombinant Human Interferon-related developmental regulator 1 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon-related developmental regulator 1 (Nerve growth factor-inducible protein PC4), IFRD1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling; Cell biology
Target Names	IFRD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.