



RPIA Antibody?FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA020104LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P49247 |
| Immunogen | Recombinant Human Ribose-5-phosphate isomerase protein (1-311AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Antigen Affinity purified & Affinity purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ribose-5-phosphate isomerase (EC 5.3.1.6) (Phosphoriboisomerase), RPIA, RPI |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | RPIA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |