



RPIA Antibody? Biotin conjugated

Product Code	CSB-PA020104LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49247
Immunogen	Recombinant Human Ribose-5-phosphate isomerase protein (1-311AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribose-5-phosphate isomerase (EC 5.3.1.6) (Phosphoriboisomerase), RPIA, RPI
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPIA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.