



RPIA Antibody? Biotin conjugated

Product Code	CSB-PA020104LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49247
Immunogen	Recombinant Human Ribose-5-phosphate isomerase protein (1-311AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribose-5-phosphate isomerase (EC 5.3.1.6) (Phosphoriboisomerase), RPIA, RPI
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPIA