



NEDD8 Antibody? Biotin conjugated

Product Code	CSB-PA618018DD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15843
Immunogen	Recombinant Human NEDD8 protein (81AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	NEDD8 (Neddylin) (Neural precursor cell expressed developmentally down-regulated protein 8) (NEDD-8) (Ubiquitin-like protein Nedd8), NEDD8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NEDD8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.