



PPARA Antibody? Biotin conjugated

Product Code	CSB-PA018421LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07869
Immunogen	Recombinant Human Peroxisome proliferator-activated receptor alpha protein (1-237AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisome proliferator-activated receptor alpha (PPAR-alpha) (Nuclear receptor subfamily 1 group C member 1), PPARA, NR1C1 PPAR
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PPARA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.