



CXCL12 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA343703LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P48061 |
| Immunogen | Recombinant Human Stromal cell-derived factor 1 protein (22-89AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Antigen Affinity purified & Affinity purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Stromal cell-derived factor 1 (SDF-1) (hSDF-1) (C-X-C motif chemokine 12) (Intercrine reduced in hepatomas) (IRH) (hIRH) (Pre-B cell growth-stimulating factor) (PBSF) [Cleaved into: SDF-1-beta(3-72), SDF-1-alpha(3-67)], CXCL12, SDF1 SDF1A SDF1B |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | CXCL12 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |