



CXCL12 Antibody? Biotin conjugated

Product Code	CSB-PA343703LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48061
Immunogen	Recombinant Human Stromal cell-derived factor 1 protein (22-89AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stromal cell-derived factor 1 (SDF-1) (hSDF-1) (C-X-C motif chemokine 12) (Intercrine reduced in hepatomas) (IRH) (hIRH) (Pre-B cell growth-stimulating factor) (PBSF) [Cleaved into: SDF-1-beta(3-72), SDF-1-alpha(3-67)], CXCL12, SDF1 SDF1A SDF1B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CXCL12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.