



FABP4 Antibody, HRP conjugated

Product Code	CSB-PA007945LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15090
Immunogen	Recombinant Human Fatty acid-binding protein adipocyte protein (6-132AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fatty acid-binding protein, adipocyte (Adipocyte lipid-binding protein) (ALBP) (Adipocyte-type fatty acid-binding protein) (A-FABP) (AFABP) (Fatty acid-binding protein 4), FABP4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FABP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.