



BAG3 Antibody, HRP conjugated

Product Code	CSB-PA002531LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95817
Immunogen	Recombinant Human BAG family molecular chaperone regulator 3 protein (220-420AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	BAG family molecular chaperone regulator 3 (BAG-3) (Bcl-2-associated athanogene 3) (Bcl-2-binding protein Bis) (Docking protein CAIR-1), BAG3, BIS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BAG3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.