



FGGY Antibody, HRP conjugated

Product Code	CSB-PA853393LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96C11
Immunogen	Recombinant Human FGGY carbohydrate kinase domain-containing protein (1-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	FGGY carbohydrate kinase domain-containing protein (EC 2.7.1.-), FGGY
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FGGY
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.