



HTR3E Antibody? Biotin conjugated

Product Code	CSB-PA010894LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A5X5Y0
Immunogen	Recombinant Human 5-hydroxytryptamine receptor 3E protein (26-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	5-hydroxytryptamine receptor 3E (5-HT3-E) (5-HT3E) (Serotonin receptor 3E), HTR3E
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HTR3E
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.