



CBY1 Antibody? Biotin conjugated

Product Code	CSB-PA896889LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3M2
Immunogen	Recombinant Human Protein chibby homolog 1 protein (41-126AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein chibby homolog 1 (ARPP-binding protein) (Cytosolic leucine-rich protein) (PIGEA-14) (PKD2 interactor, Golgi and endoplasmic reticulum-associated 1), CBY1, ARB1 C22orf2 CBY PGEA1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CBY1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.