



WISP2 Antibody?Biotin conjugated

Product Code	CSB-PA026120LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76076
Immunogen	Recombinant Human WNT1-inducible-signaling pathway protein 2 protein (24-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	WNT1-inducible-signaling pathway protein 2 (WISP-2) (CCN family member 5) (Connective tissue growth factor-like protein) (CTGF-L) (Connective tissue growth factor-related protein 58), WISP2, CCN5 CT58 CTGFL
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	WISP2