



WISP2 Antibody? Biotin conjugated

Product Code	CSB-PA026120LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76076
Immunogen	Recombinant Human WNT1-inducible-signaling pathway protein 2 protein (24-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	WNT1-inducible-signaling pathway protein 2 (WISP-2) (CCN family member 5) (Connective tissue growth factor-like protein) (CTGF-L) (Connective tissue growth factor-related protein 58), WISP2, CCN5 CT58 CTGFL
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	WISP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.