



# SATB2 Antibody? Biotin conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA892170LD01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                           |
| <b>Uniprot No.</b>         | Q9UPW6  |
| <b>Immunogen</b>           | Recombinant Human DNA-binding protein SATB2 protein (228-369AA)                         |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human   |
| <b>Tested Applications</b> | ELISA   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Biotin  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4        |
| <b>Purification Method</b> | >95%, Antigen Affinity purified & Affinity purified                                     |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | DNA-binding protein SATB2 (Special AT-rich sequence-binding protein 2), SATB2, KIAA1034 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Research Area</b>       | Signal Transduction   |
| <b>Target Names</b>        | SATB2   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.             |