



# SATB2 Antibody?Biotin conjugated

Product Code	CSB-PA892170LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPW6
Immunogen	Recombinant Human DNA-binding protein SATB2 protein (228-369AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein SATB2 (Special AT-rich sequence-binding protein 2), SATB2, KIAA1034
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SATB2