



DROSHA Antibody, HRP conjugated

Product Code	CSB-PA889096LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRR4
Immunogen	Recombinant Human Ribonuclease 3 protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Antigen Affinity purified & Affinity purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease 3 (EC 3.1.26.3) (Protein Drosha) (Ribonuclease III) (RNase III) (p241), DROSHA, RN3 RNASE3L RNASEN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DROSHA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.